

Product Information

MEM Nonessential Amino Acids (100x)

Catalogue Number: GBMNA01/01F

General Information

To encourage cell growth and survival, non-essential amino acids are added to the medium for cell culture. The recipe has 100 times more non-essential amino acids than the standard Minimum Essential Medium (MEM).

Applications:

- o As an extra source of amino acids for cells that are capable of generating non-essential amino acids
- o For cells that are capable of generating robust growth

Specification

Appearance : Clear solution

Storage : Store at +2°C to +8°C, protected from light.

Shipping Conditions: Ambient

Formulation

Components	Concentration mg/L
Glycine	0.75
L-Alanine	0.89
L-Asparagine Monohydrate	1.50
L-Aspartic Acid	1.33
L-Glutamic Acid	1.47
L-Proline	1.15
L-Serine	1.05

This product is for research use only.

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com