

## Product Information

### MEM Vitamin Solution (100X)

Catalogue Number: GBMVS01/01F

## General Information

MEM Vitamin Solution is a growth supplement that is added to cell culture medium to support cell viability and development. The Minimum Essential Medium (MEM) Vitamin Solution is made to have 100X of vitamins (MEM).

## Specification

Appearance : Clear pale yellow liquid

Storage and shelf life: Store at +2°C to +25°C upon arrival. Once opened store at +2°C to +25°C and use within 24 months

Shipping Conditions : Ambient

## Formulation

Components	Concentration mg/L
<b>Vitamins</b>	
Choline chloride	100.0
D-Calcium pantothenate	100.0
Folic Acid	100.0
Nicotinamide	100.0
Pyridoxal hydrochloride	100.0
Riboflavin	10.0
Thiamine hydrochloride	100.0
i-Inosito	200.0
<b>Inorganic Salts</b>	
Sodium Chloride (NaCl)	8500.0

**This product is for research use only.**

## Need help?

If you have any further queries, please feel free to email our cell culture specialists at [info@genexisbiotech.com](mailto:info@genexisbiotech.com)