

Product Information

Nuclease-Free Water (DEPC Treated)

Catalogue Number: GBNFWD01/01F

General Information

Nuclease-free Deionized water that has undergone DEPC treatment has been filtered into the final bottle and autoclaved. It doesn't need to be prepared, mixed, or autoclaved before use. The endonuclease, exonuclease, and RNase activity of this nuclease-free water is thoroughly tested for contamination.

Specification

Appearance : Clear liquid

Storage and shelf life : Store at 4°C.

Shipping Conditions : Ambient

This product is for research use only.

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com