

Product Information

Tris acetate-SDS Running Buffer (10X)

Catalogue Number: GBPTASB01

General Information

The Tris Acetate-SDS Running Buffer (10X) is a reliable buffer for electrophoresis of proteins using SDS-PAGE. It ensures consistent separation by maintaining an ideal pH and ionic strength in the gel system. The 10X concentrated format provides flexibility and convenience, allowing for easy dilution to a 1X working solution as needed.

Product Specification

Appearance : Clear, colorless solution

Storage & Shelf Life : Store at Room temperature and use within 12 months.

Shipping Conditions : Ambient

pH : 8.2

Concentration : 10X

Formulation

Components	Concentration g/ml
Tris base	30.29
Glacial Acetic Acid	14.43 mL
SDS	5

For Research Use Only

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.