

Product Information

RPMI 1640, with Stable Glutamine

Catalogue Number: GBRP06 / 06F

Product Specification

Appearance	: Red Colored, Clear solution
pH	: 7.00 -7.60
Storage & Shelf Life :	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping Conditions:	Ambient

Formulation

Components	Concentration mg/L	Components	Concentration mg/L
Vitamins:		Amino Acids:	
p-Amino Benzoic Acid	1.00	Glycine	10.00
D-Biotin	0.20	L-Alanyl-L-Glutamine	434.40
Choline chloride	3.00	L-Arginine	200.00
D-Calcium Pantothenate	0.25	L-Asparagine H ₂ O	56.82
Folic Acid	1.00	L-Aspartic Acid	20.00
myo-Inositol	35.00	L-Cystine 2 HCl	65.20
Nicotinamide	1.00	L-Glutamic Acid	20.00
Pyridoxal HCl	1.00	L-Histidine HCl H ₂ O	20.27
Riboflavin	0.20	L-Hydroxy-L-Proline	20.00
Thiamine HCl	1.00	L-Isoleucine	50.00
Vitamin B12	0.005	L-Leucine	50.00
Inorganic Salts:		L-Lysine HCl	40.00
Ca(NO ₃) ₂ 4 H ₂ O	100.00	L-Methionine	15.00
KCl	400.00	L-Phenylalanine	15.00
MgSO ₄ 7H ₂ O	100.00	L-Proline	20.00
NaCl	6000.00	L-Serine	30.00
NaHCO ₃	2000.00	L-Threonine	20.00
Na ₂ HPO ₄	800.00	L-Tryptophan	5.00
Other Components:		L-Tyrosine 2 Na 2H ₂ O	28.83
D-Glucose	2000.00	L-Valine	20.00
L-Glutathione Reduced	1.00		
Phenol Red Sodium Salt	5.30		

This product is for research use only.

Need help?

If you have any further queries, please feel free to email our cell culture specialists at info@genexisbiotech.com

Disclaimer: User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.